



ADVAMRYL


Poliglecaprone - 25 (Undyed & Violet)

Monofilament Synthetic Absorbable Surgical Suture

**Supreme handling performance for confidence
in subcuticular skin closure**

























**ADVA TRAY™ - the confidence to
suture with ease**



Predictable strength for superior support

Product Description




















Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
1	ADVA1642	UNDYED	0	70	1/2 CIRCLE		30	-	12
2	ADVA1742	VIOLET	0	70	1/2 CIRCLE		30	-	12
3	ADVA3442	VIOLET	0	70	1/2 CIRCLE		31	-	12
4	ADVA1006	UNDYED	1	90	1/2 CIRCLE		40	-	12
5	ADVA3759U	UNDYED	1	70	1/2 CIRCLE		40	-	12
6	ADVA1005	UNDYED	2-0	70	3/8 CIRCLE		30	-	12
7	ADVA1602	UNDYED	2-0	70	5/8 CIRCLE		26	-	12
8	ADVA1641	UNDYED	2-0	70	1/2 CIRCLE		30	-	12
9	ADVA1665	UNDYED	2-0	70	1/2 CIRCLE		26	-	12
10	ADVA1666	UNDYED	2-0	70	1/2 CIRCLE		36	-	12
11	ADVA1666S	UNDYED	2-0	45	1/2 CIRCLE		36	-	12
12	ADVA1741	VIOLET	2-0	70	1/2 CIRCLE		30	-	12
13	ADVA3441	VIOLET	2-0	70	1/2 CIRCLE		31	-	12
14	ADVA945P	UNDYED	2-0	90	1/2 CIRCLE		36	-	12
15	ADVA1001	UNDYED	3-0	70	1/2 CIRCLE		20	-	12
16	ADVA1004	UNDYED	3-0	70	3/8 CIRCLE		24	-	12
17	ADVA1008	UNDYED	3-0	90	3/8 CIRCLE		19	-	12
18	ADVA1010	UNDYED	3-0	75	3/8 CIRCLE		26	-	12
19	ADVA1326	UNDYED	3-0	70	3/8 CIRCLE		25	PRIME	12
20	ADVA1326S	UNDYED	3-0	45	3/8 CIRCLE		25	PRIME	12
21	ADVA1326MR	VIOLET	3-0	70	3/8 CIRCLE		25	PRIME	12
22	ADVA1663	UNDYED	3-0	70	1/2 CIRCLE		22	-	12

 Round Bodied

 Oval Round Bodied

 Reverse Cutting

 Cutting

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
23	ADVA1664	UNDYED	3-0	70	1/2 CIRCLE		26	-	12
24	ADVA1737	VIOLET	3-0	70	1/2 CIRCLE		20	-	12
25	ADVA3326M	UNDYED	3-0	70	3/8 CIRCLE		26	-	12
26	ADVA3326MS	VIOLET	3-0	45	3/8 CIRCLE		25	-	12
27	ADVA3665	VIOLET	3-0	75	1/2 CIRCLE		31	-	12
28	ADVA1003	UNDYED	4-0	70	3/8 CIRCLE		24	-	12
29	ADVA1007	UNDYED	4-0	75	1/2 CIRCLE		17	-	12
30	ADVA1009	UNDYED	4-0	45	3/8 CIRCLE		19	-	12
31	ADVA1205	UNDYED	4-0	70	3/8 CIRCLE		16	PRIME	12
32	ADVA1205MR	VIOLET	4-0	70	3/8 CIRCLE		16	PRIME	12
33	ADVA1648	UNDYED	4-0	70	3/8 CIRCLE		16	LONG TAPER	12
34	ADVA1748	VIOLET	4-0	70	3/8 CIRCLE		16	-	12
35	ADVA3205	UNDYED	4-0	45	3/8 CIRCLE		16	PRIME	12
36	ADVA845G	UNDYED	4-0	45	3/8 CIRCLE		16	PRIME	12
37	ADVA1002	UNDYED	5-0	45	3/8 CIRCLE		19	PRIME	12
38	ADVA3203	UNDYED	5-0	45	3/8 CIRCLE		13	PRIME	12
39	ADVA3204	UNDYED	5-0	45	3/8 CIRCLE		16	PRIME	12
40	ADVA493G	UNDYED	5-0	45	3/8 CIRCLE		13	PRIME	12
41	ADVA844G	UNDYED	5-0	45	3/8 CIRCLE		16	PRIME	12

 Round Bodied

 Oval Round Bodied

 Reverse Cutting

 Cutting

Product Description


Composition	Glycolide and ε-Caprolactone
Coating	None
Sizes	U.S.P. size 1 to 5-0
Color	Undyed and Violet
Absorption Profile	Essentially complete up to 90 days
Needle Availability	A wide variety of Round Bodied, Cutting, Reverse cutting, and Oval Round Bodied needles

Tensile Strength

Tensile Strength Retention		Days of Implantation
DYED	UNDYED	
60%	55%	7 Days (1 to 5-0)
40%	20%	14 Days (1 to 3-0 & 5-0)
40%	15%	14 Days (4-0)


66% more knot pull tensile strength Vs. In-house limit ensures excellent support throughout the critical wound healing period

Knot pull tensile strength (KgF)¹




REFERENCE _____
1 : Data on file with AMS - Bench test done on Poliglecaprone - 25, size 3-0

Monofilament Material



Violet



Undyed



Provides excellent support for a range of critical wound healing period to secure closure in a variety of tissues, procedures and patients.

AMS
Advanced MedTech Solutions

Ceejay House, 305, Level 3, Shiv Sagar Estate, Dr. Annie Besant Road, Worli, Mumbai 400018 India

T : (91)22-24946100 / 01 | amsltd.com | info@amsltd.com